



Metal -Free Digital Thermometer



KD-1370

Specifications

Model	KD-1370
Measuring Range	32.0°C ~ 43.9°C(90.0°F~109.9°F)
Resolution	0.1°C (0.1°F)
Accuracy	±0.1°C (34.0°C to 42.0°C); ±0.2°C(other range) ±0.2°F(93.2°F to 107.6°F); ±0.4°F(other range)
Display	Liquid Crystal Display 3 1/2 digits
Operation Condition	Temperature: 10 to 40°C (50 to 104°F) Relative Humidity: 15 to 95 %
Storage Condition	Temperature: -20 to 60°C (-4 to 140°F) Relative Humidity: 15 to 95 %
Battery	Micro Alkaline 192 or LR41 (1.55V)
Power Consumption	0.15 milliwatt per measurement
Low Battery Indication	Replace the battery if the low battery indicator appears
Auto Off	Automatically after 9 minutes of non use
Size	129mm x 21mm x 13mm
Weight	14g w/ battery
Standard	Complies with ASTM-E1112, EN12470-3, EN60601-1

Specifications are subject to change for improvement without notice.

Accessories:

Storage Case Probe Cover (optional)

K-jump Health Co., Ltd.

No.56 Wu Kung 5th Rd., WuKu Industrial Park,
 Taipei Hsien, Taiwan 248
 Tel: 886-2-22991378 Fax:886-2-22991386
 E-mail: service@kjump.com.tw
 http://www.kjump.com.tw



Features:



Regular



Beeper



Memory



Waterproof


Optional Functions:



Fever Alarm



C/F Switchable



HealthCare at Home
 We Care for Your Health